

Ultrasonic Beauty System

Skin Artist




Woman become vivid with fresh and elastic skin

"Skin is mirror of mind"

When Skin-cell breathe, Spoiled tissues are recovered, fine wrinkles get smoothed, sloping skin get back elastic, discoloration, dark spots, freckles, pimple will be gone.

" Now you are a skin beauty "
Ultrasonic Skin Artist which is being exported
to Japan makes you a clean, elastic beauty



▲ Recovered healthy Skin-cel tissues



▲ Ultrasonic Skin Artist




▲ Senilitied skin-cell tissues

Sloping Skin, Fine wrinkles, Fatness Diet effect.

Ultrasonic vibration with 1.7million times/sec infiltrates deeply untill inner Skin, promote to recover elastin, kolagen, so that fine wrinkles get Smoothed, get back elastic. also it promotes to resolve, ignite fat, so that changes into beautiful body.

Discoloration, dark spots, freckles pimple.

Metabolism become active by massage process, thermo process, cleaning process, so that restrain in filtration of melanin pigment which is formed by ultraviolet rays and promote to eliminate it. These processes change your skin into white, beautiful.



Find again lost yourself
- Skin Artist -

Massage process

Micro massage, which is delivered until cell unit, warms up muscles and makes your skin active. Also ultra speed massage (1.7million times/sec) infiltrates 3~4cm of innerskin and accelerates resolution of fat, so that helps your diet.



Thermo process

Ultrasonic waves heat a retired spot of your skin. then accelerates a stream of blood and lymph. Becoming normal of metabolism of cell, becomes normal and your skin will be recovered.



Cleaning process

Ultrasonic waves vibration has a function to make lots of bubble.

They remove waste matte of pores, so that makes your skin clean, soft.



Keep humidified process

Skin become rough by declining of keep humidified ability. When you use ultrawaves skin artist, it goes up, so that your skin is maintained moist and shiny.



Total Skin Care Soutlion "Skin Beau DUO"




●Can control 4 levels of output - Low, Mid, Hihg, Auto

- ◇Auto Mode automatically prepares Low, Mid, High, level of ultrasonic output by for optimum skin treatment
- ◇Ion Massager (providing 4 levels of output - Low, Mid, High, Auto)
- ◇Cleansing Massager (providing 4 levels of output - Low, Mid, High, Auto)
- ◇Vitamin Supplier (pro-vitamin mode setting)
- ◇Multiple uses of ION (Cleansing) and ultrasonic possible providing remarkable results (adopted new technology)
 - ◇Provides 4 types of Program Modes equal to these of busines uses
 - Sensitive skin treatment (Program 1)
 - Special acne treatment (Program 2)
 - Winkle Care (Program 3)
 - Freckles treatment as well as a face massage (Program 4)
 - ◇High-tech individual deauty equipment suitable for all types of skin
 - ◇Complying with maximum safety standard (safety guaranteed)
 - ◇Various skin treatment functions
 - ◇Ultrasonic beauty equipment

Ultrasonic

Ultrasonic turns into the energy when applied to inside the body, General effects of ultrasonic can be summarized as followings.

- thermal effect
- mechanical effect
- chemical effect



1) Thermal Effect

- * Ultrasonic can be absorbed into skin tissues,generating heat in skin,fat,muscles
Heat is generated in proportion to the intensity of ultasonic.
Ultrasonic thermal effect is available only in the applied skin parts.
- * Stimulating blood circulation,lymph circulation,strengthening blood vessels,boosting metabolism

2) Mechanical Effect

Ultrasonic has various mechanical effect on body parts

- * Pulse of ultrasonic gives massaging effect on skin tissues
- * Promoting the metabolism of collagen
- * Activating fibroblast (connective tissue forming cells)
- * Enhancing skin elasticity,strengthening muscles

- * Accelerating the nutrition supply to skin tissues
- * Purifying the skin cells and tissues
- * Providing elasticity to connective tissues



3) Chemical Effect

Chemical effect of ultrasonic can be summarized as followings

- * Activating skin by inducing chemical effect on collagen layers
- * Regulating pH of skin
- * Increasing the activities of APT
- * Improving the skin condition by renewing old damaged collagen layers

Ionizing Effect

- **Self-examing effect :** Supplying nutrition,maintaining the balance
- **Neutralizing effect :** Accelerating the excretion of waste,suppressing the unpleasant smell
- **Resonance effect :** Decomposition of various types of nutrition,maintaining balance accelerating metabolism
- **Moisturizing effect :** upplying nutrition,maintaining the balance



●Vitamin

Vitamin C has a negative Ion property when dissolved in water. skin has small affinity for water, thus making it hard to penetrate into the inside skin, By applying the electric power to skin exterior and using the repulsive power of ion pole, essential elements are provided to inside the skin.

●Cleansing

Cleansing function is opposite to the ionizing function. By applying the positive pole to the head, accumulated wastes between the skin tissues are excreted to the outside of skin by absorbing power of the positive pole, thus providing the skin enhancement.